

Date: Feb. 28, 2025

## **Product Technical Sheet**

### **Product Description**

Product: Anti-contamination gown, ties on neck, ties on waist, knitted cuffs

Item No.: 410914/410915


Colour: Yellow

Material: Polypropylene+ Polyethylene

Thickness: 42g PP+PE

Size: XL/2XL Please see below table

Dimension of finished product (cm)			
	410914 XL	410915 2XL	T(cm)
A	135	140	±2
B	155	160	±2
C	58	60	±1
D	28	29	±1
E	200	200	±5
F	35	35	±3



Packing: 10pcs/bag, 100pcs/ctn

Carton meas. : 55 x 35x 46 cm

### **Suggested Sterilization**

ETO: (pre-treatment cycle 16-72 hrs) (gas contact – 3 hrs/850 mg) (evaporation 48 hrs at 40°C)

### **Suggested Storage:**

Temperature : between -5°C and

40°C Humidity : below 80%

Sunlight : not to expose to direct

sunlight Shelf-Life : 5 years

**Intended Use:** It is in liquid repellent polyethylene plastic to protect clothing from splashes during medical procedures and cleaning. It can be used in hospitals, clinics, pharmaceutical, laboratories, food processing facilities and painting applications.

### **Applicable Directive/Regulation:**

Medical Device Regulation (EU) 2017/745

PPE Regulation (EU) 2016/425

### **Applicable Standards:**

ISO 9001

ISO 13485:2016

EN 14126

### **Risk Analysis**

The raw material is approved to be safe to human health (refer to the manufacturer's safety data sheet). As for the product, no hazard is foreseeable to the safety and health of user.

### **Instruction for Use:**

#### **Usage :**

1. Take it out from the package (box or polybag), then unfold it
2. Put arms into the sleeves, tie up the band at the back of neck and waist.

#### **Cautions:**

1. This product is for single use. Do not re-use it.
2. This product is for protection of part of body, not for overall protection.
3. This product is for hygienic protection, not for any other personal protection against potential risk.
4. This product is non-sterile. It must be sterilized before being used as sterile wear.
5. Do not wear it directly without any clothes. It may cause skin uncomfortably when touching skin directly.
6. Do not wear it in the place where electrostatic prevention is requested. It is not anti-static treated.
7. Choose right size according to the user's stature and figure.