

Crown Name(WH) United Co., Ltd.

Rm. 1805 Changhang Building, No. 69 Yanjiang Avenue.,

Wuhan, 430021 China

Tel: 86-27-82280908-8818 Email: sales3@crownname.com

Date: Mar. 14, 2022

# **Product Technical Sheet**


#### **Product Description**

Product: Earmuffs protection

Item No.: 356165 Colour: White Material: PP+PE Thickness: 40g

Size: Please see below table

Dimension of finished product (cm)		
	cm	Tolerance
А	26	±1
В	60	±1
С	18	±1
D	8	±1
Е	17	±1



Packing: 100pcs/bag, 1000pcs/ctn Carton meas. : 56 x 43 x 38 cm

## **Suggested Sterilization**

ETO: (pre-treatment cycle 16-72 hrs) (gas contact – 3 hrs/850 mg) (evaporation 48 hrs at 40°C)

## **Suggested Storage:**

Temperature: between -5°C and 40°C

Humidity: below 80%

Sunlight : not to expose to direct sunlight

Shelf-Life : 3 years

**Intended Use:** Used to improve the hygiene of earmuffs.

# **Applicable Directive/Regulation:**

Regulation (EU) 2016/425



Crown Name(WH) United Co., Ltd. Rm. 1805 Changhang Building, No. 69 Yanjiang Avenue., Wuhan, 430021 China

Tel: 86-27-82280908-8818 Email: sales3@crownname.com

## **Applicable Standards:**

ISO 9001:2015

ISO 10993-5:2009

ISO 10993-10:2010

# **Risk Analysis**

The raw material is approved to be safe to human health (refer to the manufacturer's safety data sheet). As for the product, no hazard is foreseeable to the safety and health of user.

#### **Instruction for Use:**

### Usage:

- 1. Take it out from the package (box or polybag), then unfold it
- 2. Open it and put on the earmuff.

## **Cautions:**

- 1. This product is for single use. Do not re-use it.
- 2. This product is for protection of part of body, not for overall protection.
- 3. This product is for hygienic protection, not for any other personal protection against potential risk.
- 4. This product is non-sterile. It must be sterilized before being used as sterile wear.
- 5. Do not wear it directly without any clothes. It may cause skin uncomfortably when touching skin directly.
- 6. Do not wear it in the place where electrostatic prevention is requested. It is not anti-static treated.
- 7. Choose right size according to the user's stature and figure.