

Date: Sep. 16, 2021

Product Technical Sheet

Product Description

Product: soft PP pants L open ankles

Item No.: 412472/412458/412455/412459/412460


Colour: blue

Material: Polypropylene

Thickness: 47g

Size: Please see below table:

Dimension of finished products(cm)						
	412472	412458	412455	412459	412460	
	S	M	L	XL	2XL	Tolerance
A	108	110	113	114	115	±2
B	47	52	58	63	68	±2
C	15	15	15	15	15	±1
D	25	25	27	27	27	±2
E	21	21	21	21	21	±2
F	77	78	79	81	83.5	±2
G	148	156	164	172	180	±2



Packing: 10pcs/bag, 50pcs/ctn

Carton dimension :

Size of carton(CM)	
S	39.5*26*27
M-L	39.5*26*28
XL-2XL	39.5*26*32

Suggested Sterilization

ETO: (pre-treatment cycle 16-72 hrs) (gas contact – 3 hrs/850 mg) (evaporation 48 hrs at 40°C)



Crown Name(WH) United Co., Ltd.

19F Bund House 7 Building, No.8 Sanyang Road, Wuhan 430010, China

Tel: 86-27-82280908-8818 Email: sales3@crowname.com

Suggested Storage:

Temperature : between -5°C and 40°C

Humidity : below 80%

Sunlight : not to expose to direct sunlight

Shelf-Life : 3 years

Intended Use: It is polypropylene material and standard fabric for non-toxic workplace. It can be used in hospital and pharmaceutical applications.

Applicable Directive/Regulation:

Directive 93/42/EEC for medical device Class I

EN 13795:2011+A1:2013

Applicable Standards:

ISO9001

ISO13485

EN13795

ISO10993

Risk Analysis

The raw material is approved to be safe to human health (refer to the manufacturer's safety data sheet).

As for the product, no hazard is foreseeable to the safety and health of user.

Instruction for Use:

Usage :

1. Take it out from the package (box or polybag), then unfold it
2. Put pants on, tie the laces on the waist.

Cautions:

1. This product is for single use. Do not re-use it.
2. This product is for protection of part of body, not for overall protection.
3. This product is for hygienic protection, not for any other personal protection against potential risk.
4. This product is non-sterile. It must be sterilized before being used as sterile wear.
5. Choose right size according to the user's stature and figure.