

Date: Mar. 14, 2022

## Product Technical Sheet

### Product Description

Product: soft PP pants L open ankles

Item No.: 412472/412458/412455/412459/412460


Colour: blue

Material: PP

Thickness: 47g

Size: Please see below table:

Dimension of finished products(cm)						
	412472	412458	412455	412459	412460	
	S	M	L	XL	2XL	Tolerance
A	108	110	113	114	115	±2
B	47	52	58	63	68	±2
C	15	15	15	15	15	±1
D	25	25	27	27	27	±2
E	21	21	21	21	21	±2
F	77	78	79	81	83.5	±2
G	148	156	164	172	180	±2



Packing: 10pcs/bag, 50pcs/ctn

Carton dimension :

Size of carton(CM)	
S	39.5*26*28
M-L	39.5*26*32
XL-2XL	39.5*26*34



Crown Name(WH) United Co., Ltd.  
Rm. 1805 Changhang Building, No. 69 Yanjiang Avenue.,  
Wuhan, 430021 China  
Tel: 86-27-82280908-8818 Email: sales3@crowname.com

---

### **Suggested Sterilization**

ETO: (pre-treatment cycle 16-72 hrs)  
(gas contact – 3 hrs/850 mg)  
(evaporation 48 hrs at 40°C)

### **Suggested Storage:**

Temperature : between -5°C and 40°C  
Humidity : below 80%  
Sunlight : not to expose to direct sunlight  
Shelf-Life : 5 years

**Intended Use:** Used to minimize the spread of infective agents to patients' surgical sites and equipment, through prevention of dispersal of bacteria-carrying skin scales from the operating room staff, thereby helping to prevent post- operative surgical site infections.

### **Applicable Directive/Regulation:**

Medical Device Regulation (EU) 2017/745

### **Applicable Standards:**

ISO13485:2016  
ISO10993-5:2019  
ISO10993-10:2019

### **Risk Analysis**

The raw material is approved to be safe to human health (refer to the manufacturer's safety data sheet).  
As for the product, no hazard is foreseeable to the safety and health of user.

### **Instruction for Use:**

#### **Usage :**

1. Take it out from the package (box or polybag), then unfold it
2. Put pants on, tie the laces on the waist.

#### **Cautions:**

1. This product is for single use. Do not re-use it.
2. This product is for protection of part of body, not for overall protection.
3. This product is for hygienic protection, not for any other personal protection against potential risk.
4. This product is non-sterile. It must be sterilized before being used as sterile wear.
5. Choose right size according to the user's stature and figure.