

Crown Name(WH) United Co., Ltd. Bund House 7 Building 19F, No.8 Sanyang Road, Wuhan 430010, China

Tel: 86-27-82280908-8818 Email: sales3@crownname.com

Date: May 24, 2024

Product Technical Sheet


Product Description

Product: SMS jacket, short sleeves

Item No.: 481677/481678/481679/481680/481681/481682/481683

Colour: green
Material: SMS
Thickness: 45gsm

Size: Please see below table



Dimension of finished product (cm)								
	481677	481678	481679	481680	481681	481682	481683	T(cm)
	XS	S	М	L	XL	2XL	3XL	1 (CIII)
Α	76	77	78	79	79	81	83	±2
В	57	60	63	66	70	74	78	±2
С	25	27	27	27	29	29	30	±2
D	24	26	27	29	29	30	30	±2
Е	19	20	21	22	23	24	25	±2

Packing: 10pcs/bag, 50pcs/ctn Carton dimension: 50*30*35cm

Suggested Sterilization

ETO: (pre-treatment cycle 16-72 hrs) (gas contact – 3 hrs/850 mg) (evaporation 48 hrs at 40°C)

Suggested Storage:

Temperature: between -5°C and 40°C

Humidity : below 80%

Sunlight : not to expose to direct sunlight

Shelf-Life : 5 years



Crown Name(WH) United Co., Ltd. Bund House 7 Building 19F, No.8 Sanyang Road, Wuhan 430010, China

Tel: 86-27-82280908-8818 Email: sales3@crownname.com

Intended Use: It is nonwoven material and standard fabric for non-toxic workplace. It can be used in hospital, clinical and pharmaceutical applications.

Applicable Directive/Regulation:

MDR for medical device Class I EN 13795:2011+A1:2013

Applicable Standards:

ISO9001 ISO13485 EN13795

Risk Analysis

The raw material is approved to be safe to human health (refer to the manufacturer's safety data sheet). As for the product, no hazard is foreseeable to the safety and health of user.

Instruction for Use:

Usage:

- 1. Take it out from the package (box or polybag), then unfold it
- 2. Put arms into the sleeves, fasten the buttons in front.

Cautions:

- 1. This product is for single use. Do not re-use it.
- 2. This product is for protection of part of body, not for overall protection.
- 3. This product is for hygienic protection, not for any other personal protection against potential risk.
- 4. This product is non-sterile. It must be sterilized before being used as sterile wear.
- 5. Choose right size according to the user's stature and figure.