



Analyticon Biotechnologies GmbH, Am Muehlenberg 10, 35104 Lichtenfels, Germany

TO WHOM IT MAY CONCERN

15th October, 2024

Declaration storage and transport temperature

We, **Analyticon Biotechnologies GmbH**, who is an established and well-known manufacturer and producer of *in vitro* Diagnostics based at Am Muehlenberg 10, 35104 Lichtenfels, Germany, hereby declare the following:

The shipping and storage conditions are the following:

Product Name	Article no.	Storage conditions	Transport conditions	Transport stability
CombiScreen® GAK	93107	2-30 °C	2-30 °C	product is stable for 3 days at 50°C / 5 days at -20°C
CombiScreen® 11SYS Plus	94100	2-30 °C	2-30 °C	product is stable for 3 days at 50°C / 5 days at -20°C
CombiScreen® 5SYS Plus	94109	2-30 °C	2-30 °C	product is stable for 3 days at 50°C / 5 days at -20°C
CombiScreen® 7SYS Plus	94110A	2-30 °C	2-30 °C	product is stable for 3 days at 50°C / 5 days at -20°C
CombiScreen® 7SYS Plus	94110	2-30 °C	2-30 °C	product is stable for 3 days at 50°C / 5 days at -20°C
CombiScreen® Nitrit	94506	2-30 °C	2-30 °C	product is stable for 3 days at 50°C / 5 days at -20°C
CombiScreen® Dip Check	93010	2-8 °C	2-8 °C*	product is stable for 18 days at ambient temperature (below 37°C)
Urilyzer® 100 Pro	UL0100Pro	5-40 °C 10-85 % rel H	-25 to +60 °C Max 75 % at 30°C rel H	n.a.

* For transport of CombiScreen® Dip Check we recommend the use of cold packs.

The products are not considered as dangerous goods.

For and on behalf of Analyticon Biotechnologies GmbH

15.10.2024

X B. Meißner

Dr. Britta Meißner
Senior Manager Regulatory Affairs & Quality A...
Signiert von: Meißner Britta