

**Product Technical Sheet**

**Product Description**

Product: TYPE 5B/6B microporous coverall with hood, knitted cuffs, elastic ankles and waist


Item No.: 410821/410822/410823/410824/410825

Colour: white

Material: Microporous film laminated

Thickness: 65g

Size: Please see below table



Dimention of finished product (CM)							
	S	M	L	XL	XXL	XXXL	T(cm)
A	161	165	170	175	180	186	±2
B	116	120	137	140	148	156	±2
C	82	84	86	87	90	93	±2
D	79	80	81	83	87	88	±2

Packing: 1pcs/bag, 50pcs/ctn

Carton measurement. : 52\*32\*58cm

**Suggested Sterilization**

ETO: (pre-treatment cycle 16-72 hrs)

(gas contact – 3 hrs/850 mg)

(evaporation 48 hrs at 40°C)

**Suggested Storage:**

Temperature : between -5°C and 40°C



**Crown Name(WH) United Co., Ltd.**  
**Rm. 1805 Changhang Building, No. 69 Yanjiang Avenue.,**  
**Wuhan, 430021 China**  
**Tel: 86-27-82280908-8818 Email: sales3@crowname.com**

---

Humidity : below 80%  
Sunlight : not to expose to direct sunlight  
Shelf-Life : 5 years

**Intended Use:** It is polypropylene material and standard fabric for non-toxic workplace. It can be used in hospitals, clinics, pharmaceutical, laboratories, food processing facilities and painting applications.

**Applicable Directive/Regulation:**

Regulation (EU) 2016/425

**Applicable Standards:**

TYPE 5B/6B

EN ISO 13688:2013

EN ISO 13982-1:2004+A1:2010

EN 13034:2005+A1:2009

EN 1073-2:2002

EN 14126:2003+AC : 2004

EN 1149-5: 2018

EN 14325:2004

**Risk Analysis**

The raw material is approved to be safe to human health (refer to the manufacturer's safety data sheet). As for the product, no hazard is foreseeable to the safety and health of user.

**Instruction for Use:**

**Usage :**

1. Take it out from the package (box or polybag), then unfold it
2. Put legs into the coverall pants, put arms into the sleeves, close the zipper and put hoods on head.

**Cautions:**

1. This product is for single use. Do not re-use it.
2. This product is for protection of part of body, not for overall protection.
3. This product is for hygienic protection, not for any other personal protection against potential risk.
4. This product is non-sterile. It must be sterilized before being used as sterile wear.
5. Choose right size according to the user's stature and figure.