



- 1** Choose proper size or cut into right size after wound debridement.



- 2** Apply dressing directly onto the wound and remove the protector film.



- 3** Put a secondary absorbent dressing on top and secure them if needed.



Type	REF Code	Size (cm)	Size (in)
Double-Sided Adhesive	SC070050DW	7 x 5	2.8" x 2"
	SC100075DW	10 x 7.5	4" x 3"
	SC160080DW	16 x 8	6.4" x 3.2"
	SC100100DW	10 x 10	4" x 4"
	SC150120DW	15 x 12	6" x 4.8"
	SC175100DW	17.5 x 10	7" x 4"
	SC150150DW	15 x 15	6" x 6"
	SC320200DW	32 x 20	12.8" x 8"
Single-Sided Adhesive	SC070050S	7 x 5	2.8" x 2"
	SC100075S	10 x 7.5	4" x 3"
	SC160080S	16 x 8	6.4" x 3.2"
	SC100100S	10 x 10	4" x 4"
	SC150120S	15 x 12	6" x 4.8"
	SC175100S	17.5 x 10	7" x 4"
	SC150150S	15 x 15	6" x 6"
	SC320200S	32 x 20	12.8" x 8"


Remark: Special specifications according to the order of the contract.

**FORYOU MEDICAL**  
**HUIZHOU FORYOU MEDICAL DEVICES CO., LTD**

HUIZHOU FORYOU MEDICAL DEVICES CO., LTD.  
ADD: No.1, Shangxia North Road, Dongjiang Hi-tech Industrial Park, Huizhou, Guangdong, China  
TEL: +86-752-5302185 FAX: +86-752-5302088 E-MAIL: medical2@foryougroup.com  
en.foryoumedical.com

LUOFUCON®

# Silicone Wound Contact Dressing



## Painless

Gentle comfort with  
painless removal.


**FORYOU MEDICAL**  
for you, for us, for life

**FORYOU MEDICAL**  
for you, for us, for life





# LUOFUCON<sup>®</sup> Silicone Wound Contact Dressing



- Red Wound
- Yellow Wound
- Pink Wound
- S: Support & Service

## Primary gentle care makes every dressing no harm at changes or removal.

LUOFUCON<sup>®</sup> Silicone Wound Contact Dressing is an open mesh, soft silicone wound contact layer that allows exudate to pass through to the secondary dressing. And is designed for the management of dry to heavily exuding wounds.



## Flexible & Soft

Able to be cut into various shapes or sizes to fit in different wound sites.



## 100% Open Passage

Open mesh enables moderate-to-high volumes of exudate to pass through into the secondary dressing.

### Protective Backing

Holds the dressing in position during application  
Ensures the dressing not sticking to gloves / hands

### Silicone Wound Contact Dressing

100% open mesh  
Distributes wound exudate fast & evenly  
Gentle and secure with pain free removal  
Without stripping of epidermal cells

### Releasing Liners

Easy peel-off liners

### Indications

LUOFUCON<sup>®</sup> Silicone Wound Contact Dressing is indicated for a wide range of wounds such as:

- Pressure ulcers
- Leg and foot ulcers
- Traumatic wounds
- Surgical wounds
- Donor sites
- Skin tears
- Burns ( I ° and II ° )

The dressing can also be used in combination with topical treatments / hydrogels and may be used in conjunction with negative pressure wound therapy.

### Contraindications

Avoid usage if ever allergic to silicone rubber.

### Benefits

- Sufficient tack for the dressing to remain in place;
- Conform to the wound and retain its position;
- Prevent pain at dressing change or removal;
- Allow free passage of exudate through the mesh;
- Comfortable & allergy-free.

