



# Katrin Hand Sanitizing Gel 1 000 ml

**Product number:** 89346

**Family:** Hand Sanitiser (Biocide)

**Number of portions:** 500 portions

**Scent/Fragrance:** not scented

**How to use:**

For hygienic hand disinfection according to EN 1500:

Highly efficient and quick-acting hand sanitising gel which spreads quickly on the skin and leaves the skin feeling clean. Mild, moisturizing and skin-friendly.



1  
Wash your hands with warm water and soap.



2  
Dry your hands with a paper towel.




3  
Use at least 3-5 ml or two clicks of the Katrin dispenser.



4  
Rub your hands for 30 seconds including your fingertips, nails and thumbs.

**Where to use:**

Hygienic Hand Disinfection in washrooms after hand washing, e.g. offices, kitchen, work places, schools, kindergartens. Suitable for food processing companies. Conforms with the European standard EN 1500.



**Product features**

**Your advantage**

Spectrum of activity ————— Bactericidal (including MRSA), Tuberculocidal, Yeasticidal, Virucidal including polio-, adeno- and noroviruses.

**Antimicrobial properties**

Hand disinfection			15 sec.	30 sec.	1 min.	1,5 min.
application recommendation for hygienic hand disinfection*	EN 1500/VAH			•		
application recommendation for surgical hand disinfection*	EN 12791 /VAH					•
tuberculocidal (M. terrae)	EN 14348	clean conditions	•			
Viruses			15 sec.	30 sec.	1 min.	1,5 min.
active against all enveloped and non-enveloped viruses (virucidal)	EN 14476	incl. HBV/HIV/HCV			•	
active against polioviruses	EN 14476	clean conditions			•	
active against adenoviruses	EN 14476	clean conditions			•	
active against noroviruses (MNV)	EN 14476	clean conditions	•			
Additional EN test results			15 sec.	30 sec.	1 min.	1,5 min.
bactericidal	EN 13727	dirty conditions	•			
yeasticidal (Candida albicans)	EN 13624	dirty conditions	•			

\* including phase 2 step 1 and phase 2 step 2 tests (quantitative suspension tests and quantitative carrier tests)

- No wastage ————— Our sealed cartridge makes it possible to use the soap right to the end. Some residual soap (less than 2% of the original volume) may be caught in the folds of the empty cartridge.
- Hygienic ————— Pump and cartridge are sealed through sonic welding preventing bacteria cross contamination.
- Efficient ————— Highly efficient and quick-acting due to more than 80 % alcohol.
- Caring ————— Contains Glycerin, which leaves the skin well-moistened and soft. Dermatologically tested and approved.



# Katrin Hand Sanitizing Gel 1 000 ml

## Product specifications

- Regulatory compliance — The biocide products commercialised under the KATRIN brand are in compliance with the Biocidal Product Regulation (EU) No 528/2012 and are approved for human hygiene purposes. The products are registered and authorized in the respective countries, where the products are sold. The product is in conformity with the REACH regulation EC/1907/2006 and its amendments. Listed in the disinfectant list of the VAH and in the virucidal product list of the IHO.
- Composition — 100 g solution contain: 80.8 g (85 % V/V) Ethanol.
- Shelf life — Unopened: please check expiry date on the label. Ones opened, use within 12 months.
- Storage — Keep between 5 – 25 degrees C during storage. Avoid direct sunlight. Keep dry.
- Handling information — Use disinfectants/ biocides safely. Always read the label and product information before use.  
Highly flammable liquid and vapour. If medical advice is needed, have product container or label at hand. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Keep container tightly closed. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. For more information, please refer to the Safety Data Sheet.

- Certificate —



ISO 9001, ISO 14001

## Suitable dispenser



**90229**  
Katrin Soap Dispenser  
White, 1 000 ml




**92209**  
Katrin Soap Dispenser  
Black, 1 000 ml

## Shipping data

### Consumer pack

X = 245 mm  
Y = 91.5 mm  
Z = 68.5 mm




Weight: 1.09 kg  
EAN: 6414301038855

### Transport pack

X = 250 mm  
Y = 150 mm  
Z = 280 mm




Weight: 6.78 kg  
EAN: 6414300089346



### Pallet

X = 1 420 mm  
Y = 800 mm  
Z = 1 200 mm



Weight: 673 kg  
EAN: 6414300089353

